

Product Datasheet

F122A Antibody Blocking peptide (orb1456529)

Description	F122A Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody AP14015c was selected from the Center region of F122A. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	Q96E09
Expiration Date	12 months from date of receipt.