

Product Datasheet

AP1G2 Antibody Blocking peptide (orb1456732)

Description	AP1G2 Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody AP13809a was selected from the N-term region of AP1G2. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	O75843
Expiration Date	12 months from date of receipt.