

Product Datasheet

LRRC37B Antibody Blocking peptide (orb1457386)

Description	LRRC37B Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody AP13189b was selected from the C-term region of LRRC37B. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	Q96QE4
Expiration Date	12 months from date of receipt.