

Product Datasheet

CHST9 Antibody Blocking peptide (orb1457387)

Description	CHST9 Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody AP13188a was selected from the N-term region of CHST9. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	Q7L1S5
Expiration Date	12 months from date of receipt.