

Product Datasheet

MLL (HRX) Antibody Blocking peptide (orb1459455)

Description	MLL (HRX) Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody was selected from the Center region of human HRX . A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	Q03164
Expiration Date	12 months from date of receipt.