



Product Datasheet

VASP (Ab57) Antibody (orb14598)

Description	VASP (Ab57) Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P50552
Expiration Date	12 months from date of receipt.



Western blot analysis of VASP(Ab-157) an...



Confocal immunofluorescent analysis of H...