



Product Datasheet

4E-BP1 (Ab-46) Antibody (orb14606)


Description	4E-BP1 (Ab-46) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q13541
Expiration Date	12 months from date of receipt.



Confocal immunofluorescent analysis of H...



Western blot analysis of 4E-BP1(Ab-45) a...



Immunohistochemical analysis of formalin...