

## Product Datasheet

### TSH-Receptor, B-Chain antibody (orb1460987)

|                            |                                                    |
|----------------------------|----------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to TSH-Receptor, B-Chain |
| <b>Species/Host</b>        | Mouse                                              |
| <b>Reactivity</b>          | Human                                              |
| <b>Conjugation</b>         | Unconjugated                                       |
| <b>Tested Applications</b> | FC, IF, IHC-P                                      |
| <b>Note</b>                | For research use only                              |
| <b>Clonality</b>           | Monoclonal                                         |
| <b>Uniprot ID</b>          | <b>P16473</b>                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                    |