


Product Datasheet

Cobicistat (GS-9350) (orb146101)

Description	Cobicistat (GS-9350) is a potent and selective inhibitor of CYP3A with IC50 of 30-285 nM.
Note	For research use only
Formula	C ₄₀ H ₅₃ N ₇ O ₅ S ₂
MW	776.023
CAS Number	[1004316-88-4]
SMILES	<chem>CN(C(N[C@@H](CCN1CCOCC1)C(N[C@@H](CC2=CC=CC=C2)CC[C@@H](NC(OCC3=CN=CS3)=O)CC4=CC=CC=C4)=O)=O)CC5=CSC(C(C)=N5</chem>
Expiration Date	12 months from date of receipt.



Chemical
structure of
Cobicistat (GS-
935...