

Product Datasheet

TPEN (orb146116)

Description	TPEN is a cell permeable, high-affinity heavy metal (Zn ²⁺ >Fe ²⁺ >Mn ²⁺) chelator that exhibits low affi...
Note	For research use only
Formula	C ₂₆ H ₂₈ N ₆
MW	424.54
CAS Number	[16858-02-9]
SMILES	<chem>N(CC1=NC=CC=C1)(CC2=NC=CC=C2)CCN(CC3=NC=CC=C3)CC4=NC=CC=C4</chem>
Expiration Date	12 months from date of receipt.