

Product Datasheet

CD70 / TNFSF7 (Activated T- & B-Lymphocyte Marker) - With BSA and Azide antibody (orb1461622)

Description	Mouse monoclonal antibody to CD70 / TNFSF7 (Activated T- & B-Lymphocyte Marker) - With BSA and Azide
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC-Fr, IP
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P32970
Expiration Date	12 months from date of receipt.