



Product Datasheet

Stathmin 1 (Ab6) Antibody (orb14617)

Description	Stathmin 1 (Ab6) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P16949
Expiration Date	12 months from date of receipt.



Western blot analysis of Stathmin 1(Ab-1...



Immunohistochemical analysis of formalin...