


Product Datasheet

FGF Receptor 1 (Ab54) Antibody (orb14619)

Description	FGF Receptor 1 (Ab54) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P08581
Expiration Date	12 months from date of receipt.



Western blot analysis of FGF Receptor 1(...)