

Product Datasheet

MAP3K1 (Mitogen-Activated Protein Kinase Kinase Kinase 1) - With BSA and Azide antibody (orb1462235)

Description	Mouse monoclonal antibody to MAP3K1 (Mitogen-Activated Protein Kinase Kinase Kinase 1) - With BSA an...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q13233
Expiration Date	12 months from date of receipt.