

## Product Datasheet

### PSMA antibody (orb1463175)

|                            |                                   |
|----------------------------|-----------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to PSMA |
| <b>Species/Host</b>        | Mouse                             |
| <b>Reactivity</b>          | Human                             |
| <b>Conjugation</b>         | Unconjugated                      |
| <b>Tested Applications</b> | IHC                               |
| <b>Target</b>              | FOLH1                             |
| <b>Note</b>                | For research use only             |
| <b>Isotype</b>             | IgG1                              |
| <b>Clonality</b>           | Monoclonal                        |
| <b>MW</b>                  | 84331                             |
| <b>Uniprot ID</b>          | <a href="#">Q04609</a>            |
| <b>Entrez</b>              | <a href="#">2346</a>              |
| <b>Dilution Range</b>      | IHC: Ready-to-use                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |