

Product Datasheet

GFP antibody (orb1463263)

Description	Goat polyclonal antibody to GFP (green fluorescent protein). GFP is a protein composed of 238 amino ...
Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, IHC-P, WB
Immunogen	Purified recombinant peptide produced in E. coli.
Target	GFP
Preservatives	NaHCO ₃ 100 mM pH 9; no preservatives
Form/Appearance	Polyclonal antibody supplied as a 200 µl (2 mg/ml) aliquot in NaHCO ₃ 100 mM pH 9. This antibody is epitope-affinity purified from goat antiserum and does not contain preservatives.
Concentration	2 mg/ml
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only
Application notes	The antibody solution should be gently mixed before use.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:500-1:5,000, IF:1:50-1:1,000, IHC-P:1:50-1:1,000, IHC-F:1:50-1:1,000, IEM:1:50-1:1,000
Expiration Date	12 months from date of receipt.