

Product Datasheet

BFP antibody (orb1463285)

Description	Goat polyclonal antibody to BFP (Blue fluorescent protein). BFP is a basic (constitutively fluoresce...
Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	IEM, IF, IHC-Fr, IHC-P, WB
Immunogen	Purified recombinant fluorescent protein produced in E. coli.
Target	BFP
Preservatives	PBS, 20% glycerol and 0.05% sodium azide
Form/Appearance	Polyclonal antibody supplied as a 100 µl (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.
Concentration	3 mg/ml
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only
Application notes	The antibody solution should be gently mixed before use.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:500-1:5,000, IF:1:50-1:500, IHC-P:1:50-1:500, IHC-F:1:50-1:500, IEM:1:50-1:500
Expiration Date	12 months from date of receipt.