

Product Datasheet

Spike RBD Domain (SARS-CoV-2) antibody (orb1463312)

Description

The spike or S glycoprotein is a

molecular weight of about 150 kDa
foun...

Species/Host

Goat

Reactivity

Virus

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

Affinity purified recombinant fusion protein using the spike protein RBD domain (residues 343 to 400) and produced in E. coli.

Target

Spike RBD Domain (SARS-CoV-2)

Preservatives

PBS, 20% glycerol and 0.05% sodium azide

Form/Appearance

Polyclonal antibody supplied as a 100 µl (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.

Concentration

3 mg/ml

Storage

Store at -20°C for long-term storage.
Store at 2-8°C for up to one month.

Note

For research use only

Application notes

The antibody solution should be gently mixed before use.

Isotype

IgG

Clonality

Polyclonal

Dilution Range

WB:1:500-1:2,000

Expiration Date

12 months from date of receipt.