

Product Datasheet

FGF8 antibody (orb146461)

Description	Goat polyclonal to FGF8
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	FGF8
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Approx 25kDa band observed in lysates of cell line Jurkat (calculated MW of 24.7 according to NP_006110.1). Recommended concentration: 1-3 µg/ml.Experiment Notes: This antibody is expected to recognize all reported isoforms (NP_149355.1; NP_006110.1; NP_149354.1; NP_149353.1)
Protein Sequence	HTTEQSLRFEFLNYP
Clonality	Polyclonal
MW	23.5; 24.7; 26.5; 27.7
Entrez	2253
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat lysate u...