


Product Datasheet

ATG4C antibody (orb146462)

Description	Goat polyclonal to ATG4C
Species/Host	Goat
Reactivity	Bovine, Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ATG4C
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Approx 50kDa band observed in Human Skeletal Muscle lysates (calculated MW of 52.5kDa according to NP_116241.2). Recommended concentration: 2-6 µg/ml. Experiment Notes: Reported variants represent identical protein: NP_116241.2, NP_835739.1
Protein Sequence	EDEKKQLKRFSTEE
Clonality	Polyclonal
MW	52.5; 52.5
Entrez	84938
Dilution Range	ELISA: 1:16000, WB: 2-6 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Skeletal ...