



Product Datasheet

CD44 antibody (orb146473)

Description	Goat polyclonal to CD44
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	CD44
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Approx 75-82kDa band observed in lysates of cell line U-2 OS (calculated MW of 81.5kDa according to NP_000601.3 and 76.6 according to NP_001001389.1). Recommended concentration: 0.01-0.03 µg/ml. Experiment Notes: This antibody is expected to recognize reported isoforms 1, 2, 3, 4, 6, 7, 8 (NP_000601.3; NP_001001389.1; NP_001001390.1; NP_001001391.1; NP_001189484.1; NP_001189485.1; NP_001189486.1).
Protein Sequence	NRDGTRYVQKGEYRT
Clonality	Polyclonal
MW	81.5
Entrez	960
Dilution Range	ELISA: 1:128000, WB: 0.01-0.03 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of U-2 OS lysate u...



Immunohistochemical staining of Human Sp...