



Product Datasheet

Cortactin (Ab-466) Antibody (orb14649)

Description	Cortactin (Ab-466) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q14247
Expiration Date	12 months from date of receipt.



Western blot analysis of Cortactin(Ab-46...



Immunohistochemical analysis of formalin...