


Product Datasheet

Mouse CCL24/Eotaxin-2 ELISA Kit (orb146692)

Description	Mouse CCL24/Eotaxin-2 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V27-V119
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	12871 MW
Uniprot ID	Q9JKC0
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
CCL24/Eotaxin-
2 PicoKine
ELISA Kit...