


Product Datasheet

Mouse CD14 ELISA Kit (orb146694)

Description	Mouse CD14 ELISA Kit
Reactivity	Mouse
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A18-P345
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	39204 MW
Uniprot ID	P10810
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 3-8 ng/ml.
Expiration Date	Please enquire.



Mouse CD14
PicoKine
ELISA Kit
standard c...