


Product Datasheet

Human TNFRSF18/GITR ELISA Kit (orb146695)

Description	Human TNFRSF18/GITR ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q26-E161
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	26000 MW
Uniprot ID	Q9Y5U5
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
TNFRSF18/GITR
PicoKine ELISA
Kit s...