


Product Datasheet

Rat RAGE ELISA Kit (orb146698)

Description	Rat RAGE ELISA Kit
Reactivity	Rat
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q24-A342
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	42664 MW
Uniprot ID	Q63495
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 78 pg/ml.
Expiration Date	Please enquire.



Rat RAGE
PicoKine
ELISA Kit
standard
cur...