


Product Datasheet

Human Progranulin ELISA Kit (orb146699)

Description	Human Progranulin ELISA Kit
Reactivity	Human
Range	1.56 ng/ml - 100 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: T18-L593
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva, urine and human milk.
MW	63544 MW
Uniprot ID	P28799
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is 33.8-44.7 ng/ml.
Expiration Date	Please enquire.



Human
Progranulin
PicoKine
ELISA Kit
sta...