


Product Datasheet

Mouse Chitinase 3-like 1/YKL-40 ELISA Kit (orb146701)

Description	Mouse Chitinase 3-like 1/YKL-40 ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Y22-A381
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	43893 MW
Uniprot ID	Q61362
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is around 100 ng/ml.
Expiration Date	Please enquire.



Mouse
Chitinase 3-
like 1/YKL-
40
PicoKine...