

Product Datasheet

Mouse TIE2/TEK ELISA Kit (orb146702)

Description	Mouse TIE2/TEK ELISA Kit
Reactivity	Mouse
Range	125 pg/ml - 8,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A23-K744
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
MW	125701 MW
Uniprot ID	Q02858
Sensitivity	10 pg/ml
Dilution Range	For serum, dilution ratio is 1:10 and concentration is around 10 ng/ml. For plasma, dilution ratio is 1:1 and concentration is 1180-3122 pg/ml.
Expiration Date	Please enquire.



Mouse
TEK/TIE2
PicoKine
ELISA Kit
standa...