


Product Datasheet

Human A2M/Alpha 2-Macroglobulin ELISA Kit (orb146708)

Description	Human A2M/Alpha 2-Macroglobulin ELISA Kit
Reactivity	Human
Range	625 pg/ml - 40,000 pg/ml
Tested Applications	ELISA
Immunogen	Target protein isolated from human plasma
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	163291 MW
Uniprot ID	P01023
Sensitivity	20 pg/ml
Dilution Range	Dilution ratio of 1:200000, concentration in serum and plasma is 0.63-1.7mg/ml.
Expiration Date	Please enquire.



Human
A2M/alpha2-
Macroglobulin
PicoKine ...