

Product Datasheet

CD307a Antibody (PE) (orb1467378)

Description

CD307a is a type I transmembrane glycoprotein of the

Reactivity

Human

Conjugation

PE

Tested

FC

Applications
Immunogen

DNA-immunization followed by a boost with the CD307a transfected cells

Target

CD307a

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.

Isotype

Mouse IgG1 kappa

Clonality

Monoclonal

Purity

Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

Q96LA6

Hazard

Non-Toxic

Information
Entrez


115350

Dilution Range


Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.

Expiration Date


12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human CD19
positive B
cell...



Flow
cytometry
surface
staining
pattern ...