

Product Datasheet

HSP70/HSC70 antibody (orb146756)

Description

Mouse monoclonal to Hsp70 (PerCP). Hsp70

	kDa hsps (hsp70s). ...
Species/Host	Mouse
Reactivity	Bovine, Canine, Drosophila, Fish, Frog, Gallus, Guinea pig, Hamster, Human, Monkey, Mouse, Other, Plant, Porcine, Rabbit, Rat, Sheep
Conjugation	PerCP
Tested Applications	FC, ICC, IEM, IHC, IP, WB
Immunogen	Recombinant HSP70/HSC70
Target	HSP70/HSC70
Preservatives	PBS pH 7.2, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SMC-104 was sufficient for detection of HSP70/HSC70 in 20 µg of heat shocked HeLa cell lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
MW	73kDa
Uniprot ID	P0DMV8 , P0DMV9
NCBI	NP_005336.3
Entrez	3303
Dilution Range	WB (1:1000), IHC (1:100), ICC/IF (1:50)
Expiration Date	12 months from date of receipt.