

Product Datasheet

GRP94 antibody [Out of stock] (orb146760)

Description

Rat monoclonal to GRP94 (ATTO-488).

	gp96) is a constitutively e...
Species/Host	Rat
Reactivity	Bovine, Canine, Equine, Frog, Gallus, Guinea pig, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat, Sheep
Conjugation	ATTO-488
Tested Applications	FC, ICC, IP, WB
Immunogen	Purified Grp94 isolated from chicken oviducts
Target	GRP94
Preservatives	PBS pH 7.2, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	0.5 µg/ml of SMC-105 was sufficient for detection of Grp94 in 20 µg of heat shocked HeLa cell lysate by colorimetric immunoblot analysis using Goat anti-rat IgG:HRP as the secondary antibody.
Isotype	IgG2a
Clonality	Monoclonal
MW	98kDa
Uniprot ID	P14625
NCBI	NP_003290.1
Entrez	7184
Dilution Range	WB (1:2000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.