

Product Datasheet

HSP90 antibody [Out of stock] (orb146793)

Description

Mouse monoclonal to Hsp90 (ATTO-390).

expressed heat...

Species/Host

Mouse

Reactivity

Canine, Fish, Gallus, Hamster, Human, Mouse, Rabbit, Rat

Conjugation

ATTO-390

Tested

AM, ELISA, ICC, IHC, IP, WB

Applications
Immunogen

Recombinant human HSP90beta

Target

HSP90

Preservatives

PBS pH 7.2, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-107 was sufficient for detection of HSP90beta in 20 µg of heat shocked HeLa cell lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG2a

Clonality

Monoclonal

MW

90kDa

Uniprot ID

[P08238](#)

NCBI

[NP_031381.2](#)

Entrez

[3326](#)

Dilution Range

WB (1:2500), IHC (1:100)

Expiration Date

12 months from date of receipt.