

Product Datasheet

HSP90 antibody (orb146808)

Description

Mouse monoclonal to Hsp90 (RPE). HSP90 is an

Species/Host	Mouse
Reactivity	Canine, Fish, Gallus, Hamster, Human, Mouse, Rabbit, Rat
Conjugation	RPE
Tested Applications	AM, ELISA, ICC, IHC, IP, WB
Immunogen	Recombinant human HSP90beta
Target	HSP90
Preservatives	PBS pH 7.2, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SMC-107 was sufficient for detection of HSP90beta in 20 µg of heat shocked HeLa cell lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG2a
Clonality	Monoclonal
MW	90kDa
Uniprot ID	P08238
NCBI	NP_031381.2
Entrez	3326
Dilution Range	WB (1:2500), IHC (1:100)
Expiration Date	12 months from date of receipt.