

Product Datasheet

HSP60 antibody [Out of stock] (orb146845)

Description

Mouse monoclonal to Hsp60 (ATTO-488).

Species/Host

Mouse

Reactivity

Bovine, Canine, Drosophila, Frog, Gallus, Guinea pig, Hamster, Human, Monkey, Mouse, Other, Plant, Porcine, Rabbit, Rat, Sheep

Conjugation

ATTO-488

Tested

ELISA, FC, IHC, IP, WB

Applications
Immunogen

Recombinant human HSP60

Target

HSP60

Preservatives

PBS, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

0.05 µg/ml of SMC-110 was sufficient for detection of HSP60 in 20 µg of heat shocked HeLa cell lysate by colorimetric immunoblot analysis using goat anti-mouse IgG as the secondary antibody.

Isotype

IgG1

Clonality

Monoclonal

MW

60kDa

Uniprot ID

[P10809](#)

NCBI

[NP_002147.2](#)

Entrez

[3329](#)

Dilution Range

WB (1:20000), IHC (1:100), ICC/IF (1:100), IP (1:200)

Expiration Date

12 months from date of receipt.