


Product Datasheet

Opioid Receptor (Ab-375) Antibody (orb14686)

Description	Opioid Receptor (Ab-375) Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P42866
Expiration Date	12 months from date of receipt.



Western blot analysis of Opioid Receptor...