



Product Datasheet

SHP-2 (Ab-542) Antibody (orb14691)


Description	SHP-2 (Ab-542) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q06124
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of SHP-2(Ab-542) a...



Confocal immunofluorescent analysis of H...