

Product Datasheet

CD74 antibody [Out of stock] (orb146932)

Description

Mouse monoclonal to CD74 (ATTO-594).

membrane protein. It...

Species/Host

Mouse

Reactivity

Human, Mouse

Conjugation

ATTO-594

Tested

FACS, FC, ICC, IHC, IP, WB

Applications
Immunogen

Human CD74 invariant chain synthetic peptide

Target

CD74

Preservatives

PBS pH 7.2, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-116 was sufficient for detection of CD74 in 20 µg of PALA cell lysates by colorimetric immunolot analysis using goat anti-mouse IgG: AP as the secondary antibody.

Isotype

IgG1

Clonality

Monoclonal

MW

33-35kDa

Uniprot ID

[P04233](#)

NCBI

[NP_001020329.1](#)

Entrez

[972](#)

Dilution Range

WB (1:1000), IHC (1:100), ICC/IF (1:50)

Expiration Date

12 months from date of receipt.