

Product Datasheet

**HSP70/HSC70 (Plant) antibody [Out of stock]
(orb146980)**

Description

Mouse monoclonal to Hsp70 (ATTO-

heat-inducible 70-kDa hsps (hsp70s...

Species/Host

Mouse

Conjugation

ATTO-390

Tested

ELISA, WB

Applications
Immunogen

Purified HSP70 from Phaseolus aureus (mung bean)

Target

HSP70/HSC70

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-120 was sufficient for detection of HSP70/HSC70 in 20 µg of Phaseolus aureus (mung bean) lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG1

Clonality

Monoclonal

MW

70kDa

Uniprot ID
[A0A1S3UWN7](#)
NCBI
[XP_014510496.1](#)
Entrez
[106769403](#)
Dilution Range

WB (1:1000)

Expiration Date

12 months from date of receipt.