

Product Datasheet

PSD95 antibody [Out of stock] (orb146998)

Description

Mouse monoclonal to PSD95 (ATTO-488).

also known as Postsyn...

Species/Host

Mouse

Reactivity

Bovine, Human, Mouse, Rat

Conjugation

ATTO-488

Tested

AM, ICC, IHC, WB

Applications
Immunogen

Recombinant rat PSD-95

Target

PSD95

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml was sufficient for detection of PSD-95 on 20 µg rat brain tissue extract by ECL immunoblot analysis using Goat Anti-Mouse IgG: HRP as the secondary.

Isotype

IgG2a

Clonality

Monoclonal

MW

75kDa

Uniprot ID

P31016

NCBI

NP_062567.1

Entrez

29495

Dilution Range

WB (1:250), IHC (1:1000), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.