

Product Datasheet

Human IGFBP-1 (C-6His) Protein (orb1471714)

Description	Recombinant Human Insulin-like Growth Factor-binding Protein 1 is produced by our Mammalian expressi...
Reactivity	Human
Conjugation	Unconjugated
Target	IGFBP-1
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Tag	C-6×His Tag
Note	For research use only
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	Mammalian
MW	26.1 KDa
Uniprot ID	P08833
Expiration Date	6 months from date of receipt.