

Product Datasheet

Human FGF-19 (N-6His) Protein (orb1471725)

Description	Recombinant Human Fibroblast Growth Factor 19 is produced by our E.coli expression system and the ta...
Reactivity	Human
Conjugation	Unconjugated
Target	FGF-19
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, 1mM EDTA, pH 8.0.
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, 1mM EDTA, pH 8.0.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Tag	N-6×His Tag
Note	For research use only
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	E. Coli
MW	23.5 KDa
Uniprot ID	O95750
Expiration Date	6 months from date of receipt.