

Product Datasheet

Human IL-20RB (C-Fc) Protein (orb1471738)

Description	Recombinant Human Interleukin-20 Receptor Subunit Beta/IL-20RB is produced by our Mammalian expressi...
Reactivity	Human
Conjugation	Unconjugated
Target	IL-20RB
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Form/Appearance	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Store at ≤-70°C, stable for 6 months after receipt.Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-Fc Tag
Note	For research use only
Purity	Greater than 85% as determined by reducing SDS-PAGE.
Source	Mammalian
MW	49.6 KDa
Uniprot ID	Q6UXL0
Expiration Date	6 months from date of receipt.