

Product Datasheet

Human IL-29 (C-10His) Protein (orb1471744)

Description	Recombinant Human Interferon Lambda-1 is produced by our Mammalian expression system and the target ...
Reactivity	Human
Conjugation	Unconjugated
Target	IL-29
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Tag	C-10×His Tag
Note	For research use only
Purity	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Source	Mammalian
MW	21.4 KDa
Uniprot ID	Q8IU54
Expiration Date	6 months from date of receipt.