

Product Datasheet

Human IFN alpha6 (C-6His) Protein (orb1471769)

Description	Recombinant Human Interferon Alpha-6 is produced by our Mammalian expression system and the target g...
Reactivity	Human
Conjugation	Unconjugated
Target	IFN
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 500mM NaCl, 1mM EDTA, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 500mM NaCl, 1mM EDTA, pH 7.4.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Tag	C-6×His Tag
Note	For research use only
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	Mammalian
MW	21.1 KDa
Uniprot ID	P05013
Expiration Date	6 months from date of receipt.