

Product Datasheet

Human IRF5 (C-6His) Protein (orb1471772)

Description	Recombinant Human Interferon Regulatory Factor 5 is produced by our E.coli expression system and the...
Reactivity	Human
Conjugation	Unconjugated
Target	IRF5
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, 30% Glycerol, 1mM DTT, pH 8.0.
Form/Appearance	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, 30% Glycerol, 1mM DTT, pH 8.0.
Storage	Store at ≤-70°C, stable for 6 months after receipt.Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-6×His Tag
Note	For research use only
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	E. Coli
MW	57.1 KDa
Uniprot ID	Q13568
Expiration Date	6 months from date of receipt.