

Product Datasheet

FKBP51 antibody [Out of stock] (orb147187)

Description

Mouse monoclonal to FKBP51 (ATTO-594).

	regulation of t...
Species/Host	Mouse
Reactivity	Canine, Hamster, Human, Mouse, Rabbit, Rat
Conjugation	ATTO-594
Tested Applications	ICC, WB
Immunogen	Synthetic peptide corresponding to the residues of human FKBP51
Target	FKBP51
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:2000 dilution was sufficient for detection of FKBP51 in ~50 µg total protein using WB analysis.
Isotype	IgG1
Clonality	Monoclonal
MW	51kDa
Uniprot ID	Q13451
NCBI	NP_001139247.1
Entrez	2289
Dilution Range	WB (1:2000), ICC/IF (1:1000)
Expiration Date	12 months from date of receipt.