

Product Datasheet

FKBP51 antibody (orb147196)

Description

Mouse monoclonal to FKBP51 (HRP). Hsp90 is

the fo...

Species/Host

Mouse

Reactivity

Canine, Hamster, Human, Mouse, Rabbit, Rat

Conjugation

HRP

Tested

ICC, WB

Applications
Immunogen

Synthetic peptide corresponding to the residues of human FKBP51

Target

FKBP51

Preservatives

PBS, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

A 1:2000 dilution was sufficient for detection of FKBP51 in ~50 µg total protein using WB analysis.

Isotype

IgG1

Clonality

Monoclonal

MW

51kDa

Uniprot ID

[Q13451](#)

NCBI

[NP_001139247.1](#)

Entrez

[2289](#)

Dilution Range

WB (1:2000), ICC/IF (1:1000)

Expiration Date

12 months from date of receipt.