

## Product Datasheet

### FKBP51 antibody (orb147198)

**Description**

Mouse monoclonal to FKBP51 (PerCP). Hsp90 is

the ...

**Species/Host**

Mouse

**Reactivity**

Canine, Hamster, Human, Mouse, Rabbit, Rat

**Conjugation**

PerCP

**Tested**

ICC, WB

**Applications**
**Immunogen**

Synthetic peptide corresponding to the residues of human FKBP51

**Target**

FKBP51

**Preservatives**

PBS, 50% glycerol, 0.09% sodium azide  
\*Storage buffer may change when conjugated

**Concentration**

1 mg/ml

**Storage**

Conjugated antibodies should be stored at 4°C

**Note**

For research use only

**Application notes**

A 1:2000 dilution was sufficient for detection of FKBP51 in ~50 µg total protein using WB analysis.

**Isotype**

IgG1

**Clonality**

Monoclonal

**MW**

51kDa

**Uniprot ID**

[Q13451](#)

**NCBI**

[NP\\_001139247.1](#)

**Entrez**

[2289](#)

**Dilution Range**

WB (1:2000), ICC/IF (1:1000)

**Expiration Date**

12 months from date of receipt.